

FKM Black

Type 1 - Standard
1/2" - 3/4"Type 1 - Standard
1" and larger

Features:	Colour code: one white dot, one yellow dot. Peroxide cured
Standards:	Meets 21CFR 177.2600. Passed USP Class VI Cytotoxicity testing

Size (inch)	Part No.	Price \$
Mini		
1/2	42MP-SFY50	17.40
3/4	42MP-SFY75	17.40
Standard		
1	40MP-SFY100	22.74
1 1/2	40MP-SFY150	21.41
2	40MP-SFY200	24.05
2 1/2	40MP-SFY250	29.44
3	40MP-SFY300	34.75
4	40MP-SFY400	40.10
6	40MP-SFY600	91.89

Price Group: DS
AV03

Silicone Clear

Type 1 - Standard
1/2" - 3/4"Type 1 - Standard
1" and larger

Features:	Colour code: no dot. Peroxide cured
Standards:	Meets 21CFR 177.2600. Passed USP Class VI Cytotoxicity testing

Size (inch)	Part No.	Price \$
Mini		
1/2	42MP-XC50	5.21
3/4	42MP-XC75	5.21
Standard		
1	40MP-XC100	4.24
1 1/2	40MP-XC150	3.61
2	40MP-XC200	4.24
2 1/2	40MP-XC250	5.90
3	40MP-XC300	6.35
4	40MP-XC400	6.51

Price Group: DS
AV03

PTFE Solid White

Type 1 - Standard
1/2" - 3/4"Type 1 - Standard
1" and larger

Application:	Not recommended for use with bolted clamps
Features:	Colour code: no dot
Standards:	Meets 21CFR 177.1550. Passed USP Class VI Cytotoxicity testing

Size (inch)	Part No.	Price \$
Mini		
1/2	42MP-G50	15.71
3/4	42MP-G75	15.71
Standard		
1	40MP-G100	26.22
1 1/2	40MP-G150	26.22
2	40MP-G200	26.22
2 1/2	40MP-G250	28.86
3	40MP-G300	28.86
4	40MP-G400	52.44

Price Group: DS
AV03

PTFE is a material subject to plastic flow. When cycled through heat and then cooled down in process lines, PTFE gaskets will not return to the original configuration and it will continue to compress as the cycling continues. This plastic flow will allow leakage if clamp torques are not maintained.